



LUTING CEMENT

SECURE LUTE RMGI

THE RESIN FORTIFIED GI LUTING CEMENT

"Post Cementation Sensitivity" may impact the success of even the best of prosthesis. Besides inducing patient discomfort and adding to patient chair-time, it may sometime warrant the replacement of an otherwise perfect prosthesis.

SECURE LUTE RMGI with its **fortified resin-glass filler matrix** helps optimizes pH and minimizes tooth desiccation, thus virtually eliminating the incidence of post cementation sensitivity. Its 'Best in Class' compressive strength (156 MPa) helps reduce marginal leakage, alongside providing unmatched physio-chemical retention.



wizdent





Secure Lute RMGI's unique creamy consistency and comfortable working time makes it ideal for luting long span bridges and cementing any type of substrate Zirconia, PFM or Metal.

Fortified resin glass filler matrix helps virtually **eliminate post cementation sensitivity.**

Increased hydrophobicity helps **minimize microleakage & prevents marginal staining**

'Best in class' compressive strength (**156 MPa**) ensures **superior adhesion & durability**

Ideal for **abutments with short height, poor retentive form & for long span bridges**

Easy **mixing & quick clean up**

High fluoride release **reduces secondary caries**





Powered by Pidilite's trusted adhesive technology, **Secure Lute RMGI** ensures stronger and more durable cementation for ALL your crowns/bridges/inlays/onlays

SPECIALLY INDICATED FOR



Cementation of Short Height Abutments / Poor Retentive Form preps & Long span bridges



Inlay/Onlay Cementation

USAGE

- ⌚ Mixing Time 30 secs
- ⌚ Working Time** 2 min 30 secs (minimum)
- ⌚ Setting Time*** 4 mins (maximum)

** Working Time is at $23\pm 1^{\circ}\text{C}$, *** setting time is at $37\pm 1^{\circ}\text{C}$

ORDERING INFORMATION

🧴 15g Powder 🧴 7ml Liquid

[Click here](#) to find more on what other dentists are talking about Wizdent Products

fb.com/fineartofdentistry

www.wizdent.com