

BONDING AGENT

# WONDER BOND SE

## THE UNIVERSAL SMART DENTAL ADHESIVE


Ever noticed your bonding agent becoming thicker with time. Most bonds tend to lose their solvent overtime due to passive evaporation within the bottle. This change in viscosity not only reduces the flow and wettability of the bonding agent but also compromises the bond strengths significantly.

**WONDER BOND SE's** unique **solvent-free** formulation addresses the issue of solvent evaporation and delivers **consistent viscosity & performance till the last drop.**



**wizdent**

  
**Pidilite**



Incorporating special **dual acid primers**, this SMART universal bond ensures deeper penetration into the dental tubules, resulting in stronger and more durable bonds.

**Consistent viscosity** and **performance** from first to last drop

Unique **All adhesive No solvent** formulation




Dual acid primers ensure **stronger** and more **durable bonds**

100% adhesive offers **more restorations per bottle**

Compatible with all etching techniques- **total-etch, self-etch & selective-etch**

**Lasting bonds** – the trusted legacy of PIDILITE



## SPECIALLY INDICATED FOR

- Adhesive for all direct composite restorations
- Adhesive for cementation of indirect restoration
- Sealant for desensitization of root surfaces

## APPLICATION

Dispense  
directly on  
applicator tip

Air thin  
for 10 sec



Scrub on moist  
tooth surface for 30 sec

Light cure for  
20 sec

## ORDERING INFORMATION

 5 ml Bottle

[Click here](#) to find more on what other dentists are talking about Wizdent Products

 7506698711

 [wizdent@pidilite.com](mailto:wizdent@pidilite.com)

 [fb.com/fineartofdentistry](https://www.facebook.com/fineartofdentistry)

 [www.wizdent.com](http://www.wizdent.com)